	Document ID: MBR-0001	Version: 1.0	Effective Date: 01AUG2021	Process Stage: Drug Product	
GO	Drug Substance Name	Batch number	Drug Substance Number	Site	
	New Product 001	6001	0.5001	USO 1	

1. Virus filtration

STEP	DESCRIPTION	RECORDED PARAMETER		RESULTS	
1.1	Obtain virus filter. Record Filter lot number and	Filter Lot Number:		XP5003T121	
1.1	expiry.	Expiry:		19 AUG 22	
1.2	Inspect the filter for irregularities	Visual Inspection		d Pass □ Fail	
1.3	Remove Product from freezer	Thaw Time Start	Date	10JUN21	
			Time	10 : 33 (AM/PM	
			Date	12JUN 21	
		Thaw Time End	Time	09: 15 (AN)/PM	
		Calculate Total Hold Days(< 7 days at Freezer)		2	
		Thaw Time within the	specification:	V Pass ⊕ Eail (1)	
1.4	Measure Conductivity	Recorded Conductivity (mS/cm)		5.5	
		Conductivity within the specification: (3.5-6.5 mS/cm)		d Pass □ Fail	
		Measure pH during the six times every hour after thaw took place. Record information on the right. (7.0-7.4)	Parameter	Value	Measurement
1.5	Measure pH		pH Measurement 1	7.3	Pass - Fail
			pH Measurement 2	7.1	Pass - Fail
			pH Measurement 3	7.0	Pass - Fail
			pH Measurement 4	7.1	Pass - Fail
			pH Measurement 5	7.1	Pass - Fail
			pH Measurement 6	7.2	Pass - Fail
		JMZ	12JUN21		
		RR	2JUNZI		

() Incorrect Entry Hold time within Range 12 JUN21